

Soil & Rock Testing Apparatus Seiken, Inc.

Specimen Consolidation Apparatus

Model No. DAC-193



It is used for consolidation of Sand specimen to do Ring shearing test.

Specification:

Specimen Ring: Initial size Dia.150X80mm, Stainless steel made

Compression method: By compressed air from 0 to 300kPa



SEIKEN, INC.

1-28-5, Higashi Ueno, Taito-ku, Tokyo 110-0015

Japan

TEL. +81-3-3834-4201 FAX.+81-3-3834-0825

URL <http://www.seikensha.com> E-mail: sks@seikensha.com